


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 11/13/17  
Time Processed : 14:19:45  
Injection Volume : 0.2  
Sample Name : Endomorphin-2  
Cat. No. : RP10926  
Lot. No. : PE4911711  
Theoretical MW : 571.67  
Observed MW : 571.45

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH