

Sample Name : Myelin Oligodendrocyte Glycoprotein Peptide (35-55), mouse, rat
 Cat. No. : RP10245
 Lot. No. : PE7301704

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump B Conc.	15	
25.00	Pumps	Pump B Conc.	75	
25.01	Pumps	Pump B Conc.	95	
30.00	Pumps	Pump B Conc.	95	
32.01	Pumps	Pump B Conc.	5	
40.00	Pumps	Pump B Conc.	5	
90.00	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	12.502	25461	3508	0.639
2	12.608	10193	1676	0.256
3	12.979	3792409	600727	95.254
4	13.262	9714	1935	0.244
5	13.491	143606	24206	3.607
Total		3981382	632053	100.000