


Mass Spectrum



Sample Information

Acquired by: Gary
Month-Day Processed: 04/03/19
Time Processed: 12:07:40 PM
Injection Volume: 0.2
Sample Name: Pancreatic Polypeptide, human
Cat. No.: RP10398
Lot. No.: PE7761904
Theoretical MW: 4181.70
Observed MW: 4182.0

Interface: ESI
Nebulizing Gas Flow: 1.5L/min
CDL Temp: 250
Block Temp: 200

Equipment: GK11010007
Interface Bias: +4.5 kV
Drying Gas Flow: 5 L/min
T.Flow: 0.2 ml/min
B.conc: 50% H2O/50% MeOH