


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 04/13/19  
Time Processed : 02:24:24 PM  
Injection Volume : 0.3  
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human  
Cat No. : RP11110  
Lot No. : PE3321904  
Theoretical MW : 2197.60  
Observed MW : 2197.5

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH