


# Mass Spectrum



## Sample Information

Acquired by : Gary  
 Month-Day Processed : 08/29/19  
 Time Processed : 06:22:01 PM  
 Injection Volume : 0.3  
 Sample Name : Myelin Basic Protein (MBP) (68-82), guinea pig  
 Cat. No. : RP10605  
 Lot. No. : E4090419H  
 Theoretical MW : 1736.80  
 Observed MW : 1736.7

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250  
 Block Temp : 200

Equipment : GK11010007  
 Interface Bias : +4.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50% H<sub>2</sub>O/50% MeOH