


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/24/18
Time Processed : 8:15:04
Injection Volume : 0.2
Sample Name : LSKL Inhibitor of Thrombospondin (TSP-1)
Cat No. : RPI9990
Lot No. : PE9371803
Theoretical MW : 459.58
Observed MW : 459.45

Interface : ESI
Nebulizing Gas Flow : 1.5 L/min
CDL Temp : 250°C
Block Temp : 200°C
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH