


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 01/26/18
Time Processed : 16:09:57
Injection Volume : 0.2
Sample Name : DYKDDDDK Peptide
Cat. No. : RP10586
Lot. No. : PE9961803
Theoretical MW : 1012.98
Observed MW : 1013.10

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH