


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 10/30/18
Time Processed : 17:49:02
Injection Volume : 0.2
Sample Name : 43Gap 26, Connexin Mimetic
Sample ID : U2845DJ250-1
Theoretical MW : 1550.78
Observed MW : 1550.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH