


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/10/16
Time Processed : 13:47:28
Injection Volume : 0.2
Sample Name : β -Amyloid (1-42), human
Cat No. : RP20527
Lot No. : PE4881611
Theoretical MW : 4514.10
Observed MW : 4513.75

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH