


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 10/16/18  
Time Processed : 17:56:45  
Injection Volume : 0.2  
Sample Name : Neurokinin B  
Cat. No. : RP10517  
Lot. No. : PE2441810  
Theoretical MW : 1210.42  
Observed MW : 1211.2

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H2O/50% MeOH