


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	05/16/19	Interface Bias	: +4.5kV
Time Processed :	03:43:47 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2ml/min
Sample Name :	Beta-Amyloid (35-25)	B.conc	: 50% $\text{H}_2\text{O}$ /50% $\text{MeOH}$
Cat No. :	RP20185		
Lot No. :	PE5091905		
Theoretical MW :	1060.30		
Observed MW :	1060.2		