


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 07/22/18  
Time Processed : 18:36:55  
Injection Volume : 0.2  
Sample Name : Ghrelin, human  
Cat No. : RPI3422  
Lot No. : PE3611807  
Theoretical MW : 3370.90  
Observed MW : 3371.40

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH