


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 05/06/20
Time Processed : 03:46:47 PM
Injection Volume : 0.3
Sample Name : Myelin Basic Protein (MBP) (synthetic)
Cat No. : RP10314
Lot No. : E0970720D
Theoretical MW : 1390.50
Observed MW : 1390.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH