


# Mass Spectrum



## Sample Information

Acquired by : Gary  
 Month-Day Processed : 10/04/19  
 Time Processed : 06:35:22 PM  
 Injection Volume : 0.4  
 Sample Name : HPV16-E711-20 epitope  
 Cat. No. : RP20285  
 Lot. No. : E7059419I  
 Theoretical MW : 1210.38  
 Observed MW : 1210.2

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250  
 Block Temp : 200

Equipment : GK11010007  
 Interface Bias : +4.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50% H2O/50% MeOH