


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/16/17
Time Processed : 22:06:25
Injection Volume : 1
Sample Name : [Leu⁵-Enkephalin
Cat No. : RPI0884
Lot No. : PE8211701
Theoretical MW : 555.62
Observed MW : 555.55

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH