


Mass Spectrum



Sample Information

Acquired by : Gary  
Month-Day Processed : 12/25/17  
Time Processed : 1:24:28  
Injection Volume : 0.4  
Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), human  
Cat. No. : RP11089  
Lot. No. : PE2361712  
Theoretical MW : 3125.59  
Observed MW : 3125.60

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias  
Drying Gas Flow  
T.Flow  
B.conc

: +4.5 kV  
: 5 L/min  
: 0.2 ml/min  
: 50%H<sub>2</sub>O/50%MeOH