


Mass Spectrum



Sample Information
Month-Day Processed : 03/10/21
Time Processed : 10:08:14 AM
Injection Volume : 0.6
Sample Name : Fibronectin CSI Peptide
Cat. No. : RP10838
Lot. No. : E8084021C
Theoretical MW : 872.96
Observed MW : 872.3

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH