

Sample Name :Kemptide
 Cat No.: RP10287
 Lot No.: E5581220F
 Time Processed :3:38:01 PM
 Month-Day-Year Processed :06/25/2020

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
31.00	Pumps	Pump A B.Conc	95
31.01	Pumps	Pump A B.Conc	5
40.00	Pumps	Pump A B.Conc	5
40.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm
 Equipment:ZJ19010325

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.603	55636	7617	0.465
2	12.936	11595366	1127393	96.881
3	13.275	66268	7670	0.554
4	13.741	54971	4132	0.459
5	14.587	117947	10830	0.985
6	16.803	78447	10235	0.655
Total		11968635	1167876	100.000