


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 11/21/18
Time Processed : 12:03:08 PM
Injection Volume : 0.2
Sample Name : GHRF (1-44), human
Cat No. : RP10734
Lot No. : PE8991811
Theoretical MW : 5039.80
Observed MW : 5040.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH