


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 11/04/18
Time Processed : 17:25:08
Injection Volume : 0.2
Sample Name : OVA G4 peptide
Cat. No. : RP20398
Lot. No. : PE4731810
Theoretical MW : 906.08
Observed MW : 906.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH