


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 06/01/19  
Time Processed : 03:51:41 PM  
Injection Volume : 0.3  
Sample Name : C-Peptide (3-33), human  
Cat No. : RP11116  
Lot No. : PE2011905  
Theoretical MW : 3020.26  
Observed MW : 3020.0

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : -3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH