


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	11/25/18	Interface Bias	: +4.5kV
Time Processed :	07:13:56 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.1	T.Flow	: 0.2ml/min
Sample Name :	[Arg⁸]-Vasotocin	B.conc	: 50%H ₂ O/50%MeOH
Cat No. :	RP10061		
Lot No. :	PE6351811		
Theoretical MW :	1050.22		
Observed MW :	1050.2		