


Mass Spectrum



Sample Information

Month-Day Processed : 08/11/20
Time Processed : 20:33:54
Injection Volume : 0.3
Sample Name : Tuftsin
Cat No. : RP10624
Lot No. : E6858020H
Theoretical MW : 500.60
Observed MW : 500.6

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH