


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/30/18
Time Processed : 18:41:51
Injection Volume : 0.2
Sample Name : BNP (1-32), rat
Cat No. : RP11121
Lot No. : PE2881809
Theoretical MW : 3452.94
Observed MW : 3452.5

Interface

:ESI
Nebulizing Gas Flow :1.5L/min
CDL Temp :250°C
Block Temp :200°C

Equipment

: GK11010007
Interface Bias :+4.5kV
Drying Gas Flow :5 L/min
T.Flow :0.2ml/min
B.conc :50%H2O/50%MeOH