


# Mass Spectrum



## Sample Information

Month-Day Processed : 10/24/20  
 Time Processed : 20:11:50  
 Injection Volume : 0.3  
 Sample Name : Brain Natriuretic Peptide (BNP) (1-32), rat  
 Cat. No. : RP11121  
 Lot. No. : E2241220J  
 Theoretical MW : 3452.94  
 Observed MW : 3452.5

## Interface

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250  
 Block Temp : 200

## Equipment

Equipment : GK11010007  
 Interface Bias : +4.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50% H2O/50% MeOH