


Mass Spectrum



Sample Information
Month-Day Processed : 10/30/20
Time Processed : 02:16:11 PM
Injection Volume : 0.4
Sample Name : Brain Natriuretic Peptide (BNP) (1-32), rat
Cat. No. : RP11121
Lot. No. : E6127820J
Theoretical MW : 3452.94
Observed MW : 3452.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH