


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/27/19
Time Processed : 02:00:26 PM
Injection Volume : 0.2
Sample Name : Pancreatic Polypeptide,human
Cat No. : RP10398
Lot No. : PE7641903
Theoretical MW : 4181.70
Observed MW : 4181.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH