


Mass Spectrum



Sample Information		Interface		Equipment	
Month-Day Processed :	02/05/21	Nebulizing Gas Flow : 1.5L/min		: GK11010007	
Time Processed :	10:08:28 AM	CDL Temp : 250		: +3.5 kV	
Injection Volume :	0.4	Block Temp : 200		: 5 L/min	
Sample Name :	Secretin, canine			Drying Gas Flow	
Cat No. :	RP12943			T.Flow	
Lot No. :	E637621A			: 0.2 ml/min	
Theoretical MW :	3069.45			B.conc	
Observed MW :	3068.8			: 50%H2O/50%MeOH	