


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	10/10/19	Interface Bias	:+4.5kV
Time Processed :	19:57:25	Drying Gas Flow	:5 L/min
Injection Volume :	0.4	T.Flow	:0.2ml/min
Sample Name :	Human Papillomavirus (HPV) E7 protein (49 - 57)	B.conc	:50%H2O/50%MeOH
Cat No. :	RP20249		
Lot No. :	E5569719J		
Theoretical MW :	1120.27		
Observed MW :	1120.0		