


Mass Spectrum



Sample Information
Month-Day Processed : 09/28/20
Time Processed : 08:55:17 PM
Injection Volume : 0.3
Sample Name : HPV16 E7(86-93)
Cat. No. : RP20499
Lot. No. : E4947820I
Theoretical MW : 815.10
Observed MW : 814.3

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH