


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 12/29/18
Time Processed : 07:43:33 PM
Injection Volume : 0.2
Sample Name : Angiotensin II (3-8), human
Cat No. : RP11253
Lot No. : PE4861812
Theoretical MW : 774.90
Observed MW : 775.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH