


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 06/11/20
 Time Processed : 09:11:49 PM
 Injection Volume : 0.3
 Sample Name : PACAP (1-27), human, ovine, rat
 Cat No. : RP10334
 Lot No. : E1500420F
 Theoretical MW : 3147.62
 Observed MW : 3146.8

Interface

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH