


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 05/11/18
Time Processed : 9:41:58
Injection Volume : 0.2
Sample Name : Rabies Virus Glycoprotein
Cat. No. : RP20464
Lot. No. : PE0191805
Theoretical MW : 3266.70
Observed MW : 3266.75

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH