


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 05/23/17
Time Processed : 10:12:01
Injection Volume : 0.2
Sample Name : Kemptide
Cat No. : RP10287
Lot No. : PE6641705
Theoretical MW : 771.90
Observed MW : 771.70

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH