


Mass Spectrum



Sample Information
User : Gary
Month-Day Processed : 08/05/20
Time Processed : 09:09:43 PM
Injection Volume : 0.3
Sample Name : CLIP (86-100)
Cat No. : RP20279
Lot No. : E6889620H
Theoretical MW : 1674.10
Observed MW : 1673.8

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200