


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 09/25/19
 Time Processed : 04:59:01 PM
 Injection Volume : 0.5
 Sample Name : Tyrosine Protein Kinase JAK 2 (pTyr8, 9)
 Cat. No. : RP20119
 Lot. No. : E48301191
 Theoretical MW : 2082.11
 Observed MW : 2082.0

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H2O/50% MeOH