


Mass Spectrum



Sample Information

| | | | |
|-----------------------|--------------------------|-----------------|---------------------------------|
| User : | Gary | Equipment | : GK11010007 |
| Month-Day Processed : | 09/19/19 | Interface Bias | : +4.5kV |
| Time Processed : | 08:42:17 AM | Drying Gas Flow | : 5 L/min |
| Injection Volume : | 0.3 | T.Flow | : 0.2ml/min |
| Sample Name : | [Met5]-Enkephalin, amide | B.conc | : 50% H ₂ O/50% MeOH |
| Cat No. : | RP10887 | | |
| Lot No. : | E4771219I | | |
| Theoretical MW : | 572.68 | | |
| Observed MW : | 572.2 | | |