


MS Spectrum Graph



Sample Information

Acquired by : Admin
Date Acquired : 2017-1-5 14:17:13
Injection Volume : 1
Sample Name : PLP(139-151)
Cat No. : RP10550
Lot No. : PE8481701
Theoretical MW : 1521.74
Observed MW : 1521.20

Interface :ESI

Interface Bias : +4.5 kV
Nebulizing Gas Flow : 1.5L/min
Drying Gas Flow : 5 L/min
CDL Temp : 250°C
T.Flow : 0.2 ml/min
Block Temp : 200°C
B.conc : 50%H2O/50%MeOH