


# Mass Spectrum



## Sample Information

Acquired by: Gary  
Month-Day Processed: 04/02/19  
Time Processed: 05:29:01 PM  
Injection Volume: 0.2  
Sample Name: Melanin Concentrating Hormone, salmon  
Cat No.: RP10662  
Lot No.: PE5811903  
Theoretical MW: 2099.48  
Observed MW: 2099.4

Interface: ESI  
Nebulizing Gas Flow: 1.5L/min  
CDL Temp: 250  
Block Temp: 200

Equipment: GK11010007  
Interface Bias: +4.5 kV  
Drying Gas Flow: 5 L/min  
T.Flow: 0.2 ml/min  
B.conc: 50% H2O/50% MeOH