


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/07/20
Time Processed : 04:57:16 PM
Injection Volume : 0.3
Sample Name : Tuftsin
Cat No. : RP10624
Lot No. : E8629020A
Theoretical MW : 500.60
Observed MW : 500.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH