


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 04/08/19
Time Processed : 02:51:09 PM
Injection Volume : 0.2
Sample Name : Melanotan (MT)-II
Cat No. : RP17323
Lot No. : PE8711903
Theoretical MW : 1024.19
Observed MW : 1024.0

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200