


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 08/31/19
Time Processed : 01:38:05 PM
Injection Volume : 0.5
Sample Name : C-Peptide, dog
Cat No. : RP11115
Lot No. : E4657619H
Theoretical MW : 3174.49
Observed MW : 3174.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH