


# Mass Spectrum



Sample Information  
Month-Day Processed : 02/06/21  
Time Processed : 9:15:19 AM  
Injection Volume : 0.6  
Sample Name : LL-37, scrambled  
Cat No. : RP20332  
Lot No. : E3760821B  
Theoretical MW : 4493.30  
Observed MW : 4492.6

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH