


# Mass Spectrum



Sample Information  
Month-Day Processed : 10/06/20  
Time Processed : 12:37:28 PM  
Injection Volume : 0.4  
Sample Name : Biotin-TAT (47-57)  
Cat No. : RP20408  
Lot No. : E35365201  
Theoretical MW : 1786.20  
Observed MW : 1786.0

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200  
Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H2O/50% MeOH