

Sample Name : Influenza NP (147-155)  
 Cat No. : RP20260  
 Lot No. : E5456421A  
 Time Processed : 20:52:30  
 Month-Day-Year Processed : 01/25/2021  
 Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow: 1 ml/min  
 Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
31.00	Pumps	Pump A B.Conc	5
31.01	Pumps	Pump A B.Conc	5
40.00	Controller	Stop	


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm  
 Equipment : ZJ19010140

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	13.075	19912	3141	0.126
2	13.286	20996	2175	0.133
3	13.653	149686	17984	0.948
4	13.759	68228	13160	0.432
5	13.954	170834	34194	1.082
6	14.136	15259935	1331041	96.643
7	14.388	68595	19827	0.434
8	14.606	31901	3081	0.202
Total		15790087	1424604	100.000