


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 07/31/17  
Time Processed : 9:09:37  
Injection Volume : 0.2  
Sample Name : [Met<sup>5</sup>-Enkephalin  
Cat No. : RPI0886  
Lot No. : PE3291707  
Theoretical MW : 573.66  
Observed MW : 573.30

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH