


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 11/03/19
Time Processed : 11:14:29 AM
Injection Volume : 0.3
Sample Name : PUMA BH3
Cat. No. : RP20510
Lot. No. : E5595419
Theoretical MW : 3049.30
Observed MW : 3048.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH