


Mass Spectrum



Sample Information	Interface	Interface Bias
Acquired by : Gary	: ESI	: +4.5 kV
Month-Day Processed : 03/21/17	Nebulizing Gas Flow : 1.5L/min	Drying Gas Flow : 5 L/min
Time Processed : 17:30:26	CDL Temp : 250°C	T.Flow : 0.2 ml/min
Injection Volume : 0.2	Block Temp : 200°C	B.conc : 50%H2O/50%MeOH
Sample Name : [Met⁵-Enkephalin		
Cat. No. : RP10886		
Lot. No. : PE4181703		
Theoretical MW : 573.66		
Observed MW : 573.25		