


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/12/19
Time Processed : 09:03:11 AM
Injection Volume : 0.4
Sample Name : CLIP (86-100)
Cat No. : RP20279
Lot No. : E1587619I
Theoretical MW : 1674.10
Observed MW : 1673.8

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200