


Mass Spectrum



Sample Information
Month-Day Processed : 01/27/21
Time Processed : 3:41:11 PM
Injection Volume : 0.3
Sample Name : Palmitoyl Tetrapeptide-3
Sample ID : U6182GA220-1
Theoretical MW : 694.91
Observed MW : 694.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH